

Product Datasheet

BGLF4 Antibody - C-terminal region (orb2245295)

Description

BGLF4 Antibody - C-terminal region

Species/Host

Rabbit

Conjugation

Unconjugated

Tested

IHC

Applications
Immunogen

This BGLF4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 386-417 amino acids from the C-terminal region of human BGLF4.

Preservatives

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

Ig

Clonality

Polyclonal

MW

48351 Da

Uniprot ID
P13288
NCBI
YP_401688.1
Expiration Date

12 months from date of receipt.