


Product Datasheet

Retinoic Acid Receptor beta antibody (orb224546)

Description	Rabbit polyclonal antibody to Retinoic Acid Receptor beta
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cell lysa...