

Product Datasheet

FGFR1 Antibody (Phospho-Tyr307) (orb2245820)

Description

FGFR1 Antibody (Phospho-Tyr307)

Species/Host

Rabbit

Conjugation

Unconjugated

Tested

DOT

Applications
Immunogen

This FGFR1 Antibody is generated from rabbits immunized with a KLH conjugated synthetic phosphopeptide corresponding to amino acid residues surrounding Y307 of human FGFR1.

Preservatives

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

Application Info: DB: 1:500

Isotype

Ig

Clonality

Polyclonal

MW

91868 Da

Uniprot ID
P11362
NCBI
NP_001167534.1
Expiration Date

12 months from date of receipt.