


Product Datasheet

RANKL antibody (orb224599)

Description	Rabbit polyclonal antibody to RANKL
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:200, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HuvEc cell lysa...