


Product Datasheet

CD271 antibody (orb224643)

Description	Rabbit polyclonal antibody to CD271
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat and NIH-3T3 cells.