


Product Datasheet

TNF-R2 antibody (orb224648)

Description	Rabbit polyclonal antibody to TNF-R2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
LOVO cell
lysate...