

Product Datasheet

NHSL1 Antibody - middle region (orb2246655)

Description	NHSL1 Antibody - middle region
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This NHSL1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 632-661 amino acids from the Central region of human NHSL1.
Preservatives	Liquid PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid PBS with 0.09% (W/V) sodium azide.
Storage	Maintain refrigerated at 2-8C for up to 2 weeks. For long term storage store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Application Info: WB: 1:1000
Isotype	Ig
Clonality	Polyclonal
MW	171 kDa
Uniprot ID	Q5SYE7
NCBI	NP_001137532.1
Expiration Date	12 months from date of receipt.