



Product Datasheet

NCAM1 antibody (orb22474)

Description	Goat polyclonal antibody to NCAM1
Species/Host	Goat
Reactivity	Bovine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Target	CD56 / NCAM1 (isoform 2)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 170kDa band observed in Human Brain (Cerebellum, Cerebral and Frontal Cortex) lysates (calculated MW of 94.6kDa according to NP_851996.2). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	KASWTRPEKQETLD
Clonality	Polyclonal
MW	94.6
Entrez	4684
Dilution Range	ELISA: 1:32000, WB: 0.01-0.03 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb22474 (0.01ug/ml)
staining of Human F...



Immunohistochemical
staining of Human
Ad...