


Product Datasheet

AES antibody (orb22478)

Description	Goat Polyclonal to AES
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	AES
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 30kDa band observed in Human Brain (Cerebellum and Cerebral Cortex) lysates (calculated MW of 29.1kDa according to NP_945320.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	SHLPQQLKFTTSD
Clonality	Polyclonal
MW	29.1; 22.0; 21.8
Entrez	166
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb22478
(0.3ug/ml)
staining of
Human Ce...