



Product Datasheet

TXNDC5 antibody (orb22482)


Description	Goat Polyclonal to TXNDC5
Species/Host	Goat
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	TXNDC5 (aa281-295)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa plus faint 35kDa bands observed in Human Lymph Nodes lysates (calculated MW of 47.6kDa according to NP_110437.2) and of 36.2kDa according to NP_001139021.1. Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	RDLESLREYVESQLQ
Clonality	Polyclonal
MW	47.6; 36.2
Entrez	81567
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Lymph Node...



Western blot analysis of Mouse Lymph Node...



Western blot analysis of Mouse Lymph Node...