


Product Datasheet

PAPD5 antibody (orb22483)

Description	Goat Polyclonal to PAPD5
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PAPD5 (aa78-90)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 75-80kDa band observed in Human and Mouse Brain lysates (calculated MW of 75.8kDa according to Human NP_001035374.2 and 74.2kDa according to Mouse NP_001157969.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	HALPAEQRDFLP
Clonality	Polyclonal
MW	75.8; 70.8
Entrez	64282
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of Hum...