

Product Datasheet

ZNF623 Antibody - C-terminal region (orb2248340)

Description

ZNF623 Antibody - C-terminal region

Species/Host

Rabbit

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This ZNF623 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 506-535 amino acids from the C-terminal region of human ZNF623.

Preservatives

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

Application Info: WB: 1:1000FC: 1:10~50

Isotype

Ig

Clonality

Polyclonal

MW

61392 Da

Uniprot ID
075123
NCBI
NP_001075949.1
Expiration Date

12 months from date of receipt.