



Product Datasheet

GNB2L1 antibody (orb22484)


Description	Goat Polyclonal to GNB2L1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	GNB2L1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 33kDa band observed in Human Lymph Nodes and Mouse and Rat Spleen lysates (calculated MW of 35.1kDa according to NP_006089.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	NTLGVCKYTVQDESH
Clonality	Polyclonal
MW	35.1
Entrez	10399
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of Hum...



Western blot analysis of staining of NIH...



Immunohistochemical staining of Human Sm...