

Product Datasheet

UNQ6493 Antibody - C-terminal region (orb2248589)

Description

UNQ6493 Antibody - C-terminal region

Species/Host

Rabbit

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This UNQ6493 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 78-107 amino acids from the C-terminal region of human UNQ6493.

Preservatives

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -2°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

Application Info: WB: 1:1000

Isotype

Ig

Clonality

Polyclonal

MW

12794 Da

Uniprot ID
Q6UXR8
Expiration Date

12 months from date of receipt.