


Product Datasheet

EPHA2/3/4 antibody (orb224907)

Description	Rabbit polyclonal antibody to EPHA2/3/4
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cell lysa...