



Product Datasheet

IDH3B antibody (orb22493)


Description	Goat polyclonal antibody to IDH3B
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IDH3B (aa33-46)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 39kDa band observed in Human, Mouse and Rat Heart and Skeletal Muscle lysates (calculated MW of 41.9kDa according to NP_777280.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	HAASRSQAEDVRVE
Clonality	Polyclonal
MW	42.2; 41.9
Entrez	3420
Dilution Range	ELISA: 1:1000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Skeletal ...



Western blot analysis of staining of Pig...



Western blot analysis of staining of Mou...