


Product Datasheet

C22ORF28 antibody (orb22496)

Description	Goat Polyclonal to C22ORF28
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	HSPC117 (aa125-138)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 55kDa band observed in lysates of cell line Jurkat (calculated MW of 55.2kDa according to NP_055121.1). Recommended concentration: 0.05-0.2 µg/ml.
Protein Sequence	RLLRTNLDESDVQP
Clonality	Polyclonal
MW	55.2
Entrez	51493
Dilution Range	ELISA: 1:64000, WB: 0.05-0.2 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb22496
(0.05ug/ml)
staining of
Jurkat ...