


Product Datasheet

HRH2 antibody (orb22498)

Description	Goat Polyclonal to HRH2
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide with sequence C-SHNSHKTSRLNNS, from the internal region of the protein sequence according to NP_001399939.1; NP_001010973.1.
Target	Histamine Receptor H2 (aa309-323) (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 45kDa band observed in Human Duodenum and Ileum lysates (calculated MW of 44.5kDa according to NP_001124527.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	SHNSHKTSRLNNS
Clonality	Polyclonal
MW	44.5; 40.1
Entrez	3274
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb22498
(0.1µg/ml)
staining of
Duodenum...