



Product Datasheet

TRAP antibody (orb2250)


Description	TRAP antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human TRAP/Tartrate Resistant Acid Phosphatase (101-200/325 aa)
Target	ACP5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34 kDa
Uniprot ID	P13686
Dilution Range	WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Rat lung...



Immunohistochemical staining of mouse em...



Immunohistochemical staining of Rat lung...