


Product Datasheet

IL7 Receptor antibody (orb22503)

Description	Goat polyclonal antibody to IL7 Receptor
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	CD127 / IL7R (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band observed in lysates of cell line NIH3T3 (calculated MW of 51.6kDa according to NP_032398.3). Recommended concentration: 2-4 µg/ml.
Protein Sequence	DDHSFWCHSQLE
Clonality	Polyclonal
MW	51.6
Dilution Range	ELISA: 1:32000, WB: 2-4 µg/ml
Expiration Date	12 months from date of receipt.



orb22503
(2ug/ml)
staining of
NIH3T3 lys...