


Product Datasheet

CB1 antibody (orb22505)

Description	Goat polyclonal antibody to CB1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CB1 (isoform a)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 52kDa band observed in lysates of cell lines A549 and NIH3T3 (calculated MW of 52.9kDa according to Human NP_057167.2 and 52.8kDa according to Mouse NP_031752.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	SNDIQYEDIKGDMS
Clonality	Polyclonal
MW	52.9; 49.1
Entrez	1268
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of A54...