


Product Datasheet

MTHFD1L antibody (orb22509)


Description	Goat polyclonal antibody to MTHFD1L
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	MTHFD1L (aa535-538)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 110kDa band observed in lysates of cell line HeLa (calculated MW of 106kDa according to NP_056255.2). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	DTNDRFLRKITIGQ
Clonality	Polyclonal
MW	106
Entrez	25902
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



WB analysis of
Human Cerebellum
lysate u...



Immunohistochemical
staining of Human
To...



Western blot analysis
of NIH3T3 (A), Mou...