



Product Datasheet

ACPP antibody (orb22510)

Description	Goat Polyclonal to ACPP
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ACPP / PAP (aa269-282)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band observed in Human Prostate lysates (calculated MW of 48.3kDa according to NP_001127666.1). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	NHMKRATQIPSYKK
Clonality	Polyclonal
MW	48.3; 44.6
Entrez	55
Dilution Range	ELISA: 1:32000, WB: 0.01-0.03 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human prostate ...



Immunohistochemical analysis of formalin...