

Product Datasheet

YAP1 antibody (orb225103)

Description	Mouse monoclonal antibody to YAP1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, IHC-P, WB
Note	For research use only
Clonality	Monoclonal
Dilution Range	ELISA: 1-10000, WB: 1-500-1:2000, IHC-P: 1-200-1:1000, FACS: 1-200-1:400
Expiration Date	12 months from date of receipt.




Figure 2: Western blot analysis using YAP1 mouse mAb against Hela (1), C6 (2) and Cos7 (3) cell lysate.

Western blot analysis of recombinant pro...