


Product Datasheet

TFB2M antibody (orb22514)

Description	Goat polyclonal antibody to TFB2M
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	TFB2M
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 45kDa band observed in Human Lung lysates (calculated MW of 45.3kDa according to NP_071761.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	KHCFGRRSATVIDH
Clonality	Polyclonal
MW	45.3
Entrez	64216
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Lung lysa...