


## Product Datasheet

### Osteopontin antibody (orb225195)

|                            |                                                                             |
|----------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to Osteopontin                                    |
| <b>Species/Host</b>        | Mouse                                                                       |
| <b>Reactivity</b>          | Human                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                |
| <b>Tested Applications</b> | ELISA, ICC, IF, IHC-P, WB                                                   |
| <b>Note</b>                | For research use only                                                       |
| <b>Clonality</b>           | Monoclonal                                                                  |
| <b>Dilution Range</b>      | ELISA: 1-10000, WB: 1-500-1:2000, IHC-P: 1-200-1:1000, IF/ICC: 1-200-1:1000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                             |



Western blot analysis of recombinant pro...