



Product Datasheet

IRF2BP1 antibody (orb22523)

Description	Goat Polyclonal to IRF2BP1
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IRF2BP1 (aa109-122)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 75kDa band observed in lysates of cell line HepG2 (calculated MW of 61.7kDa according to NP_056464.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	QDRYDRATSSGRLP
Clonality	Polyclonal
MW	61.7
Entrez	26145
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysate US...



Western blot analysis of NIH3T3 lysate u...