

Product Datasheet

TNFR1 antibody (orb225252)

Description	Rabbit polyclonal antibody to TNFR1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P: 1:50-1:200, WB: 1:500-2000, IF/ICC: 1:100-1:500, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.