



Product Datasheet

GNAQ antibody (orb22529)


Description	Goat polyclonal antibody to GNAQ
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	GNAQ / ALPHA-q (aa162-175)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 38 to 40kDa band observed in Human, Mouse and Rat Testis lysates (calculated MW of 42.1kDa according to Human NP_002063.2 and to Rat NP_112298.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	NDLDRVADPAYLPT
Clonality	Polyclonal
MW	42.1
Entrez	2776
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of Per...



Western blot analysis of staining of Mou...



Western blot analysis of staining of Pig...