



Product Datasheet

SDF4 antibody (orb22530)

Description	Goat Polyclonal to SDF4
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	SDF4 (aa161-175)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms (NP_057631.1; NP_057260.2).WB: Approx 45kDa band observed in lysates of cell line Jurkat (calculated MW of 41.8kDa according to NP_057260.2). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	EYKVKFLASKGHSEK
Clonality	Polyclonal
MW	39.6; 41.8
Entrez	51150
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml, IHC-P: 5-10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate u...



Immunohistochemical staining of Human Ki...