



Product Datasheet

ZEB2 antibody (orb22539)

Description	Goat polyclonal antibody against ZEB2
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ZEB2 (aa545-558)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 100+38kDa bands observed in Human Liver lysates (calculated MW of 136kDa according to NP_055610.1). Recommended concentration: 0.2-0.5 µg/ml. We believe the two observed bands represent two breakdown fragments.
Protein Sequence	TDSRRQISNIKKEK
Clonality	Polyclonal
MW	136; 134
Entrez	9839
Dilution Range	ELISA: 1:32000, WB: 0.2-0.5 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...



Immunohistochemical staining of Human Ki...