

## Product Datasheet

### GASP1 antibody (orb225475)

|                            |                                        |
|----------------------------|----------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to GASP1    |
| <b>Species/Host</b>        | Rabbit                                 |
| <b>Reactivity</b>          | Human                                  |
| <b>Conjugation</b>         | Unconjugated                           |
| <b>Tested Applications</b> | ICC, IF, IHC-P                         |
| <b>Note</b>                | For research use only                  |
| <b>Clonality</b>           | Polyclonal                             |
| <b>Dilution Range</b>      | IHC-P: 1:50-1:200, IF/ICC: 1:100-1:500 |
| <b>Expiration Date</b>     | 12 months from date of receipt.        |