



Product Datasheet

Cyb5r3 antibody (orb22552)

Description	Goat polyclonal antibody to Cyb5r3
Species/Host	Goat
Reactivity	Bovine, Canine, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CYB5R3 / Dia 1 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 30kDa band observed in Mouse and Rat Lung lysates (calculated MW of 34.1kDa according to NP_084063.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	CLPNLERVGHPKER
Clonality	Polyclonal
MW	34.1
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Lung lysa...



Immunohistochemical staining of Human Li...