



Product Datasheet

IRF2BP1 antibody (orb22554)

Description	Goat polyclonal antibody to IRF2BP1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IRF2BP1 (aa461-474)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60kDa band observed in Human Liver lysates (calculated MW of 61.7kDa according to NP_056464.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	QPPTQHRLVARNGE
Clonality	Polyclonal
MW	61.7
Entrez	26145
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH3T3 lysate a...



Western blot analysis of Human Liver lys...