


Product Datasheet

PAX1 antibody (orb22555)

Description	Goat polyclonal antibody to PAX1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PAX1 (aa318-329)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 55.5kDa according to NP_006183.2). Recommended concentration: 0.2-0.6 µg/ml.
Protein Sequence	PKMEDWAGVNRT
Clonality	Polyclonal
MW	55.5
Entrez	5075
Dilution Range	ELISA: 1:32000, WB: 0.2-0.6 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Cerebellu...