

Product Datasheet

Acetyl-Histone H3 (Lys23) antibody (orb225550)

Description	Rabbit polyclonal antibody to Acetyl-Histone H3 (Lys23)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.