



Product Datasheet

PDLIM1 antibody (orb225557)


Description	Rabbit polyclonal antibody to PDLIM1
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rabbit, Rat, Sheep, Xenopus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IHC-P: 1:200, IF/ICC: 1:200-500, ELISA: 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells, C...



Western blot analysis of extracts from r...



Western blot analysis of extracts from 2...