

Product Datasheet

CD28 antibody (orb22561)

Description	Goat Polyclonal to CD28
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Cd28 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 40kDa band observed in Mouse Thymus lysates (calculated MW of 25.2kDa according to NP_031668.3). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	KEKHLCHTQSSPK
Clonality	Polyclonal
MW	25.2
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Thymus ly...