


Product Datasheet

ABCC4 antibody (orb22562)

Description	Goat polyclonal antibody to ABCC4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ABCC4 (aa70-82)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 100kDa band observed in Human Lung and Human Prostate lysates (calculated MW of 96.8kDa according to NP_001098985.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	RAENDAQKPSLTR
Clonality	Polyclonal
MW	150; 96.8
Entrez	10257
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Prostate ...