


Product Datasheet

TNFRSF6B antibody (orb225625)

Description	Rabbit polyclonal antibody to TNFRSF6B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using TNFRSF6B antibody (orb225625).