


Product Datasheet

DYR1A antibody (orb225823)

Description	Rabbit polyclonal antibody to DYR1A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using DYR1A antibody (orb225823).