


Product Datasheet

TNFAIP8L2 antibody (orb225879)

Description	Rabbit polyclonal antibody to TNFAIP8L2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HUVEC cells using TNFAIP8L2 antibody (orb225879).