



Product Datasheet

DNMT3B antibody (orb226011)


Description	Rabbit polyclonal antibody to DNMT3B
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-200, IF/ICC: 1:100-500, ELISA: 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis of HT-29 cells using DNMT3B antibody (orb226011).



Immunofluorescence analysis of HT29 cells using DNMT3B antibody (orb226011).



Immunohistochemical staining of Mouse tissue using DNMT3B antibody (orb226011).