


## Product Datasheet

### BAIAP2L1 antibody (orb226371)

|                            |                                        |
|----------------------------|----------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to BAIAP2L1 |
| <b>Species/Host</b>        | Rabbit                                 |
| <b>Reactivity</b>          | Human, Mouse, Rat                      |
| <b>Conjugation</b>         | Unconjugated                           |
| <b>Tested Applications</b> | WB                                     |
| <b>Note</b>                | For research use only                  |
| <b>Clonality</b>           | Polyclonal                             |
| <b>Dilution Range</b>      | WB: 1:500-1:1000                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.        |



Western blot analysis of HepG2 cells using BAIAP2L1 antibody.