


Product Datasheet

ARPP21 antibody (orb226385)

Description	Rabbit polyclonal antibody to ARPP21
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using ARPP21 antibody (orb226385).