


## Product Datasheet

### MRP9 antibody (orb226714)

|                            |                                    |
|----------------------------|------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to MRP9 |
| <b>Species/Host</b>        | Rabbit                             |
| <b>Reactivity</b>          | Human                              |
| <b>Conjugation</b>         | Unconjugated                       |
| <b>Tested Applications</b> | WB                                 |
| <b>Note</b>                | For research use only              |
| <b>Clonality</b>           | Polyclonal                         |
| <b>Dilution Range</b>      | WB: 1:500-1:1000                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.    |



Western blot analysis of HUVEC cells 293...