


Product Datasheet

NUTF2 antibody (orb226774)

Description	Rabbit polyclonal antibody to NUTF2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using NUTF2 antibody (orb226774).