



Product Datasheet

ZDHHC7 antibody (orb226790)

Description	Rabbit polyclonal antibody to ZDHHC7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1000, IF: 1:100-500, P-ELISA: 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis of HUVEC cells using ZDHHC7 antibody.



IF image of HUVEC cells stained with ZDHHC7 antibody (DF4276).