

## Product Datasheet

### IL4 antibody (orb22681)

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Description</b>       | Mouse monoclonal antibody to IL4           |
| <b>Species/Host</b>      | Rat                                        |
| <b>Reactivity</b>        | Mouse                                      |
| <b>Conjugation</b>       | Unconjugated                               |
| <b>Note</b>              | For research use only                      |
| <b>Application notes</b> | Monoclonal Antibodies, Anti-Mouse Cytokine |
| <b>Clonality</b>         | Monoclonal                                 |
| <b>Expiration Date</b>   | 12 months from date of receipt.            |