

## Product Datasheet

### RAP1GAP antibody (orb226901)

|                            |                                       |
|----------------------------|---------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to RAP1GAP |
| <b>Species/Host</b>        | Rabbit                                |
| <b>Reactivity</b>          | Human, Mouse                          |
| <b>Conjugation</b>         | Unconjugated                          |
| <b>Tested Applications</b> | IHC-P, WB                             |
| <b>Note</b>                | For research use only                 |
| <b>Clonality</b>           | Polyclonal                            |
| <b>Dilution Range</b>      | WB: 1:500-1:1000, IHC-P: 1:50-1:200   |
| <b>Expiration Date</b>     | 12 months from date of receipt.       |